



eDNAAnalyzer MANUAL

eDNAAnalyzer is a user-friendly and open-access computational tool developed to process and filter taxonomic assignment data from metabarcoding studies, particularly derived from environmental DNA (eDNA) and invertebrate-derived DNA (iDNA) approaches.

How to cite eDNAAnalyzer? Olimpio, L.W.G.F.; Gestich, C.C.; Saranholi, B.H.; Galetti Jr, P.M.; Freitas, P.D. 2025. eDNAAnalyzer: a fast and user-friendly computational tool for processing massive taxonomic assignment data derived from eDNA and iDNA metabarcoding (doi:).

How to access the eDNAAnalyzer? You can find the eDNAAnalyzer tool by accessing the GitHub repository at <https://github.com/Leo-9821/eDNAAnalyzer>. The executable file (.exe) is available for Windows®, while the Python source code (.py files) can be used to run the program on Linux® or macOS®, in this case download the “main.py” and “metabar.py” files and “img” folder to the same directory in your computer, and make sure to install pandas, Pillow and openpyxl libraries.

How does eDNAAnalyzer work? The software provides two main functions: “Threshold application” and “Results consolidation”, each requiring a specific input file. Examples of input and output files are available in the repository https://github.com/Leo-9821/eDNAAnalyzer/tree/master/example_files.

An overview of the [general operation](#) of the program, showing its interfaces, is provided at the end of this manual, along with more detailed initial instructions for performing the two main functions of eDNAAnalyzer: “[threshold application](#)” and “[results consolidation](#)”.

Choosing the option **Threshold Application**

This option processes data from prior taxonomic assignments by calculating the total number of reads per OTU per sequencing sample, filtering out OTUs according to an adopted threshold value (default ≥ 0.05), and generating outputs following selected filters.

Process Overview:

1. Calculate the total number of reads per OTU per sequencing sample, based on previous taxonomic assignments provided ([input file](#)).
2. Filter OTUs according to an adopted cutoff value (default: $\geq 0.05\%$).
3. Generate files ([output files](#)) containing the total number of reads per OTU, the OTUs that did not pass by the threshold value and the OTUs that passed.

Input File: a file in .xlsx or .csv format, (e.g., example_input.xlsx available at https://github.com/Leo-9821/eDNAalyzer/blob/master/example_files/eng/xlsx/1_threshold_application/input/input_example.xlsx), containing the taxonomic assignment data, must be provided with the following required columns:

- sequencing_sample (identification of sequencing samples)
- area_sampler (indicating the collection area and DNA sampler - e.g., water sample, soil sample, invertebrate used as sampler, etc. - and ensure that both areas and samplers are unambiguously identified to avoid errors during separation. Use an underscore to separate the area from the sampler: area_sampler)
- point (indicating the collection point)
- n_reads (number of reads per OTU)
- taxon (taxon determined from taxonomic assignment)

Extra columns containing additional information can also be included in the table, as shown in Table 1.

Table 1: Summary table of an input file showing the general information to perform the threshold application step.

| sequencing_sample | barcode | tag | area_sampler | point | aliquot | otu | n_reads | %_id | taxon |
|-------------------|---------|-----|--------------|-------|---------|-----|---------|-------|----------------------------------|
| P08 | 12SrRNA | TA | P1_MQ | 1 | 1 | 1 | 233 | 100 | <i>Nycticorax nycticorax</i> |
| P08 | 12SrRNA | TB | P2_MC | 2 | 1 | 1 | 50 | 99.26 | <i>Bos taurus</i> |
| P08 | 12SrRNA | TB | P2_MC | 2 | 1 | 2 | 41 | 99.25 | <i>Canis lupus</i> |
| P08 | 12SrRNA | TB | P2_MC | 2 | 1 | 4 | 3 | 100 | <i>Sus scrofa</i> |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 2 | 20 | 99.26 | <i>Canis lupus familiaris</i> |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 3 | 38 | 100 | <i>Canis lupus familiaris</i> |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 4 | 524 | 100 | <i>Canis lupus familiaris</i> |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 6 | 7 | 100 | <i>Canis lupus</i> |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 7 | 2 | 99.26 | <i>Equus caballus</i> |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 8 | 4 | 100 | <i>Leopardus pardalis</i> |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 9 | 2 | 100 | <i>Puma concolor</i> |
| P08 | 12SrRNA | TF | P2_MC | 2 | 1 | 1 | 6 | 100 | plant junction region |
| P08 | 12SrRNA | TG | P2_MC | 2 | 1 | 1 | 17 | 99.29 | <i>Canis lupus</i> |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 1 | 109 | 100 | <i>Bos taurus</i> |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 2 | 52 | 100 | <i>Canis lupus familiaris</i> |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 3 | 9 | 100 | <i>Pecari tajacu</i> |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 4 | 1859 | 100 | <i>Rhea americana</i> |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 5 | 7 | 100 | <i>Didelphis albiventris</i> |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 6 | 4 | 98.56 | <i>Gallus gallus</i> |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 7 | 6 | 97.81 | <i>Gymnogyps californianus</i> |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 9 | 4 | 100 | <i>Streptopelia decaocto</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 1 | 14 | 100 | <i>Bos taurus</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 2 | 127 | 100 | <i>Equus caballus</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 3 | 28 | 100 | <i>Sus scrofa</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 4 | 27 | 100 | <i>Tamandua tetradactyla</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 5 | 8 | 100 | <i>Hydrochoerus hydrochaeris</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 6 | 25 | 100 | <i>Bubulcus ibis</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 7 | 1013 | 100 | <i>Gallus gallus</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 8 | 4 | 100 | <i>Gallus gallus</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 9 | 2 | 100 | <i>Nycticorax nycticorax</i> |
| P08 | 12SrRNA | TK | P3_MC | 3 | 1 | 1 | 33 | 100 | <i>Zaedyus pichiy</i> |
| P08 | 12SrRNA | TK | P3_MC | 3 | 1 | 2 | 26 | 97.04 | <i>Euphractus sexcinctus</i> |
| P08 | 12SrRNA | TL | P3_MC | 3 | 1 | 1 | 1111 | 100 | <i>Canis aureus</i> |
| P08 | 12SrRNA | TL | P3_MC | 3 | 1 | 2 | 481 | 100 | <i>Zaedyus pichiy</i> |

Output Files: All tables can be saved in .xlsx and .csv format. See an example of a file provided in this step in Table 2.

Table 2: Summary table of an output file showing the data processed after running the threshold application step.

| sequencing_sample | barcode | tag | area_sampler | point | aliquot | otu | n_reads | %_id | taxon | final_otu_curated |
|-------------------|---------|-----|--------------|-------|---------|-----|---------|-------|-------------------------------|-------------------|
| P08 | 12SrRNA | TA | P1_MQ | 1 | 1 | 1 | 233 | 100 | <i>Nycticorax nycticorax</i> | |
| P08 | 12SrRNA | TB | P2_MC | 2 | 1 | 1 | 50 | 99.26 | <i>Bos taurus</i> | |
| P08 | 12SrRNA | TB | P2_MC | 2 | 1 | 2 | 41 | 99.25 | <i>Canis lupus</i> | |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 2 | 20 | 99.26 | <i>Canis lupus familiaris</i> | |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 3 | 38 | 100 | <i>Canis lupus familiaris</i> | |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 4 | 524 | 100 | <i>Canis lupus familiaris</i> | |
| P08 | 12SrRNA | TG | P2_MC | 2 | 1 | 1 | 17 | 99.29 | <i>Canis lupus</i> | |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 1 | 109 | 100 | <i>Bos taurus</i> | |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 2 | 52 | 100 | <i>Canis lupus familiaris</i> | |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 4 | 1859 | 100 | <i>Rhea americana</i> | |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 1 | 14 | 100 | <i>Bos taurus</i> | |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 2 | 127 | 100 | <i>Equus caballus</i> | |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 3 | 28 | 100 | <i>Sus scrofa</i> | |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 4 | 27 | 100 | <i>Tamandua tetradactyla</i> | |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 6 | 25 | 100 | <i>Bubulcus ibis</i> | |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 7 | 1013 | 100 | <i>Gallus gallus</i> | |
| P08 | 12SrRNA | TK | P3_MC | 3 | 1 | 1 | 33 | 100 | <i>Zaedyus pichiy</i> | |
| P08 | 12SrRNA | TK | P3_MC | 3 | 1 | 2 | 26 | 97.04 | <i>Euphractus sexcinctus</i> | |
| P08 | 12SrRNA | TL | P3_MC | 3 | 1 | 1 | 1111 | 100 | <i>Canis aureus</i> | |
| P08 | 12SrRNA | TL | P3_MC | 3 | 1 | 2 | 481 | 100 | <i>Zaedyus pichiy</i> | |

Choosing the option **Results Consolidation**

This process will provide tables with the results, lists of species with their number of detections and reads. The tables can be generated by separating the lists by sampler, area, or both.

Process Overview:

1. Edit the table containing the OTUs of interest by completing the `final_otu_curated` column with the selected taxonomic assignment after manual curation. This step aims to review and refine taxonomic identification to correct possible inconsistencies in assignment; for this, data on species distribution and field observations must be incorporated for better accuracy.
2. Input the curated table ([input file](#)) into eDNAAnalyzer, and the program will process the data and then return the information of interest according to the filters selected by the user (e.g.; sampler, area, or sampler and area).
3. Generate files with consolidated results ([output files](#)) in .xlsx and .csv formats. When necessary, .csv files in a .zip file will be provided.

Input File: The input file (Table 3) for the “results consolidation” step is a table (in .xlsx or .csv format) containing the final OTU data from the previous threshold application step and the curated OTU list.

Table 3. Summary table of an input file showing the general information to perform the consolidation step.

| sequencing_sample | barcode | tag | area_sampler | point | aliquot | otu | n_reads | %_id | taxon | final_otu_curated |
|-------------------|---------|-----|--------------|-------|---------|-----|---------|-------|-------------------------------|-------------------------------|
| P08 | 12SrRNA | TA | P1_MQ | 1 | 1 | 1 | 233 | 100 | <i>Nycticorax nycticorax</i> | <i>Nycticorax nycticorax</i> |
| P08 | 12SrRNA | TB | P2_MC | 2 | 1 | 1 | 50 | 99.26 | <i>Bos taurus</i> | <i>Bos taurus</i> |
| P08 | 12SrRNA | TB | P2_MC | 2 | 1 | 2 | 41 | 99.25 | <i>Canis lupus</i> | <i>Canis lupus familiaris</i> |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 2 | 20 | 99.26 | <i>Canis lupus familiaris</i> | <i>Canis lupus familiaris</i> |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 3 | 38 | 100 | <i>Canis lupus familiaris</i> | <i>Canis lupus familiaris</i> |
| P08 | 12SrRNA | TC | P2_MC | 2 | 1 | 4 | 524 | 100 | <i>Canis lupus familiaris</i> | <i>Canis lupus familiaris</i> |
| P08 | 12SrRNA | TG | P2_MC | 2 | 1 | 1 | 17 | 99.29 | <i>Canis lupus</i> | <i>Canis lupus familiaris</i> |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 1 | 109 | 100 | <i>Bos taurus</i> | <i>Bos taurus</i> |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 2 | 52 | 100 | <i>Canis lupus familiaris</i> | <i>Canis lupus familiaris</i> |
| P08 | 12SrRNA | TH | P3_MC | 3 | 1 | 4 | 1859 | 100 | <i>Rhea americana</i> | <i>Rhea americana</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 1 | 14 | 100 | <i>Bos taurus</i> | <i>Bos taurus</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 2 | 127 | 100 | <i>Equus caballus</i> | <i>Equus caballus</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 3 | 28 | 100 | <i>Sus scrofa</i> | <i>Sus scrofa</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 4 | 27 | 100 | <i>Tamandua tetradactyla</i> | <i>Tamandua tetradactyla</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 6 | 25 | 100 | <i>Bubulcus ibis</i> | <i>Bubulcus ibis</i> |
| P08 | 12SrRNA | TI | P3_MQ | 3 | 1 | 7 | 1013 | 100 | <i>Gallus gallus</i> | <i>Gallus gallus</i> |
| P08 | 12SrRNA | TK | P3_MC | 3 | 1 | 1 | 33 | 100 | <i>Zaedyus pichiy</i> | <i>Euphractus sexcinctus</i> |
| P08 | 12SrRNA | TK | P3_MC | 3 | 1 | 2 | 26 | 97.04 | <i>Euphractus sexcinctus</i> | <i>Euphractus sexcinctus</i> |
| P08 | 12SrRNA | TL | P3_MC | 3 | 1 | 1 | 1111 | 100 | <i>Canis aureus</i> | <i>Canis lupus familiaris</i> |
| P08 | 12SrRNA | TL | P3_MC | 3 | 1 | 2 | 481 | 100 | <i>Zaedyus pichiy</i> | <i>Euphractus sexcinctus</i> |

Output File: all tables can be saved in .xlsx and .csv formats. See an example of a file provided in this step in Table 4.


Table 4. Summary table of an output file showing the consolidated results after running the consolidation step. Note that the taxon list is presented according to the filter selected in different guides of the datasheet.

| | A | B | C | D |
|----|----|----------------------------------|-------|------------------|
| 1 | | Taxon | Reads | Detections in P6 |
| 2 | 0 | <i>Tapirus terrestris</i> | 4195 | 1 |
| 3 | 1 | <i>Canis lupus familiaris</i> | 1269 | 1 |
| 4 | 2 | <i>Pitheciidae</i> | 736 | 1 |
| 5 | 3 | <i>Coendou insidiosus</i> | 451 | 1 |
| 6 | 4 | <i>Callicebus nigrifrons</i> | 313 | 1 |
| 7 | 5 | <i>Bos taurus</i> | 212 | 1 |
| 8 | 6 | <i>Lycalopex vettulus</i> | 202 | 1 |
| 9 | 7 | <i>Rattus rattus</i> | 177 | 1 |
| 10 | 8 | <i>Puma concolor</i> | 164 | 1 |
| 11 | 9 | <i>Callithrix penicillata</i> | 118 | 1 |
| 12 | 10 | <i>Columba livia</i> | 113 | 1 |
| 13 | 11 | <i>Gallus gallus</i> | 98 | 1 |
| 14 | 12 | <i>Sus scrofa</i> | 94 | 1 |
| 15 | 13 | <i>Cingulata</i> | 88 | 1 |
| 16 | 14 | <i>Pecari tajacu</i> | 65 | 1 |
| 17 | 15 | <i>Thraupis sayaca</i> | 65 | 1 |
| 18 | 16 | <i>Callithrix sp.</i> | 63 | 1 |
| 19 | 17 | <i>Euphractus sexcinctus</i> | 48 | 1 |
| 20 | 18 | <i>Cairina moschata</i> | 39 | 1 |
| 21 | 19 | <i>Hydrochaeris hydrochaeris</i> | 38 | 1 |
| 22 | 20 | <i>Canidae</i> | 36 | 1 |
| 23 | 21 | <i>Dasyus novemcinctus</i> | 32 | 1 |
| 24 | 22 | <i>Chrysocyon brachyurus</i> | 19 | 1 |
| 25 | 23 | <i>Nymphicus hollandicus</i> | 18 | 1 |
| 26 | | | | |


General Operation of the eDNAAnalyzer Program

After accessing the program at <https://github.com/Leo-9821/eDNAAnalyzer>, follow the steps by clicking on the options available.


Step 1: Select the language




Step 2: Read the manual and then choose one of the two processing options




Step 3: Choose the option threshold application



Step 4: Upload a file containing the taxonomic assignment data




Step 5: Check the provided folder directory to ensure that the file was uploaded correctly.




The screenshot shows a window titled "Threshold application" with a purple header. Inside, there is a section titled "Threshold application" with a back arrow icon. Below this is a button labeled "Upload a file". Underneath the button, it says "Uploaded file C:/Users/input_example.xlsx". There is a section titled "Threshold selection" with the text "Indicate the threshold in percentage: (default value 0.05)" followed by an empty text input field and a "%" symbol. At the bottom is a large button labeled "RUN".

Step 6: Insert a threshold value in percentage or keep empty to use the default (0.05%) and run the Threshold application step.




The screenshot shows the same "Threshold application" window. The "Uploaded file C:/Users/input_example.xlsx" is still displayed. In the "Threshold selection" section, the text input field now contains the value "0.07". The "RUN" button remains at the bottom.

Step 7: After running the threshold application, three output files will be generated, each containing, respectively: (i) the calculated threshold value per sequencing sample, (ii) the OTUs with number of reads below the adopted threshold value, and (iii) the OTUs with number of reads equal to or greater than the adopted threshold. Save at least this last file in .xlsx or .csv formats, to use it for manual curation and then as an input file for data filtering in the results consolidation step.



Step 8: Choose the option results consolidation.




Step 9: Upload a file containing a curated table after running the threshold application step

The screenshot shows a web application window titled "Results consolidation". It features a navigation bar with a back arrow and the title. Below the navigation bar, there is a section for file upload with a button labeled "Upload a file" and a status message "No file uploaded". A "Data filtering" section contains a label "Filter table by:" followed by two radio buttons: "Sampler" (which is selected) and "Area". Below this is a "Samplers definition" section with a text prompt "Insert the samplers ID as identified in the table (one per row):" and an empty text input field. At the bottom of the form is a large "RUN" button.

Step 11: Check the provided folder directory to ensure that the file was uploaded correctly.


This screenshot shows the same "Results consolidation" application window, but now a file has been uploaded. The "Upload a file" button is still present, but the status message has changed to "Uploaded file C:/Users/processed_results_curated.xlsx". The "Data filtering" section remains the same with "Sampler" selected. The "Samplers definition" section and the "RUN" button are also visible and unchanged from the previous state.

Step 12: Select the options to filter the data by sampler and/or area



The screenshot shows a web application window titled "Results consolidation". It features a navigation bar with a back arrow and the title. Below the navigation bar, there is a section for file upload with a button "Upload a file" and a message "Uploaded file C:/Users/processed_results_curated.xlsx". The "Data filtering" section has a label "Filter table by:" followed by two checked checkboxes: "Sampler" and "Area". The "Samplers definition" section contains a text input field with the placeholder text "Insert the samplers ID as identified in the table (one per row):". At the bottom, there is a large "RUN" button.

Step 12: If you selected the "sampler" option, type the sampler IDs, separating them from each other, by pressing the key "enter" from your keyboard. If you did not select the "Sampler" option in the previous step, leave this field blank.



This screenshot shows the same "Results consolidation" window, but the "Samplers definition" text input field is now populated with the text "MC" and "MQ" on two separate lines. All other elements, including the file upload section and the "Data filtering" checkboxes, remain the same as in the previous screenshot.

Step 13: Save the tables generated with the data selected according to the chosen filters. To compress all the generated tables in .csv format into a single .zip file, select the ZIP option. Thus, you will save the consolidated results with the final list of taxa, and information about the number of reads that each taxon presented and the number of times it was detected.

